



The smart choice for Lipo-sculpture in esthetic surgery field.


ULTRA-Z






Ultrasonic Cavitation

The high ultrasonic energy of ULTRA-Z probes causes micro bubbles within the tumescent solution to expand and collapse, dislodging fatty tissue and mixing them into the infiltrated solution through three steps - Creation of bubble nucleus, Growth of bubble and Implosion of bubble.



Emulsification

The phenomenon in which unwanted fatty tissue is changed to a frothy liquid form such like nectar by ultrasonic energy.



Features


- Ultra-Z is designed to optimize every stage of the body and face sculpting and emulsifies unwanted fatty tissue, while protect important tissues. Quick recovery of patient thanks to minimal invasive surgery.
- Appropriate probes with various diameters to treat all areas of the body and face

 Various ranges of diameters, 1.2mm to 3.7mm, provide tailor-made lipo-sculpture for small volumes of fatty tissues 1.2mm ULTRA-Z probe, unique and thinnest in the world, enables lipo-sculpture for face with minimal swelling, bruising, bleeding and maximal safety.



The ULTRA-Z System :::





Specification

Rated Voltage	AC110V or AC 230V	Ultrasonic Frequency	26~50 KHz	
Rated Frequency	50/60Hz	Size(W*D*H)	397mm*360mm*170mm	
Power Consumption	300VA	Weight	15 Kg	
Fuse	T3.15AL 250Vac	Output Mode	N mode(continuous) 30~40 KHz	
Protection class	Class1, Type BF	Output Mode	Z mode(pulsed) 30~40 KHz	
Leakage current	In acc. With IEC60601-1	Output W	Max 90W	

Excellent assistant for Lipo-sculpture by Ultra-Z

ZP-1000

Feature

> Irrigation

- Easy to set FLOW and RPM value using UP and DOWN buttons.
- Total volume of the injected fluid is displayed.
- Easy to reset total volume injected using CLEAR button.
- ERROR indicator lights on and warning sound indicate it when error occurs.

➤ Suction

- Oilless-piston type vacuum pump that has a strong suction power is adopted.
- Triple overflow protection system:
 - Floating valve, Overflow protection sensor, Overflow protection filter
- Silicon hose for infiltration and suction
 - Almost no transformation for long-time using.
 - No blockage thanks to good elasticity of silicon hose.
- 3,000cc bottle for suction: High strength and high thermal resistance.



Specification

Dimensions(W*D*H)				
Weight				
Power Supply		AC 110V or 230V		
Frequency				
Power Consumption	Max. 780 VA			
Fuse	T 6.3 AL/250V			
Irrigation		Suction		
Flow rate	12.5~1375 ml/min.	The degree of vacuum	680 mmHg	
Accuracy	± 10% ml/min.	Accuracy	Max. >= 680 mmHg -10%	
		Displacement	100 L/min.	
		Pump	Oil-less piston typemotor pump	

ZP-1000 for Irrigation only!



Specification

Dimensions(W*D*H)	210mm*360mm*133mm		
Weight	6Kg		
Power Supply	AC 100V ~240V		
Frequency	50/60 Hz		
Power Consumption	Max. 90 VA		
Fuse	T 1.0 AL/250V		
Flow rate	12.5~1375 ml/min.		
Accuracy	± 10% ml/min.		



Picture	Code No.	Name of Item	Specification	ULTRA-Z	ZP-1000	ZP-1000 (Irrigation only)
0.=	FS01-02R	Single Foot Switch	3 pin plug, cable 4m	0		0
	HP01-01R	Handpiece	36kHz, 159 * 30mm	0		
0=	CA05-01R	Handpiece Cable	2 pin plug 3m	0		
	SP01-01R	Spanner	Al, 35*33mm	0		
all a	PB01-01R	Probe	Ti-Alloy, 2.9 * 275mm	0		
-	PB01-03R	Probe	Ti-Alloy, 1.9 * 134mm	0		
#3	PB01-04R	Probe	Ti-Alloy, 1.2 * 65mm	0		
	PB01-05R	Probe	Ti-Alloy, 3.7 * 295mm	0		
0	FS02-03R	Double Foot Switch	4 pin plug, cable 4m		0	
0	TU01-01R	Suction Tube	Ø8mm*1.5m		0	
.0	TU01-02R	IrrigationTubing Set	Ø6mm∗4m, Needle Spike, Connector		0	0
	BT01-01R	Suction Bottle	3,000 сс		0	
+	FL01-01R	Overflow Protection Filter	Filter		0	
791	DS01-01R	Overflow Protection Sensor	Infrared Detector		0	
-	IC01-01R	Infiltration Cannula	ø2mm*20cm		0	0
•	IC01-02R	Infiltration Cannula	ø2mm*30cm		0	0
-	IC01-03R	Infiltration Cannula	ø0.8mm*7.5cm		0	0
•	IC01-04R	Infiltration Cannula	ø0.9mm*7.5cm		0	0
	IH01-01R	Infiltration Handle	luer lock, Stainless Steel		0	0
	SC01-01R	Liposuction Cannula	ø3mm*30cm, blunt type		0	
•	SC01-02R	Liposuction Cannula	ø4mm*30cm, blunt type		0	
	SH01-01R	Liposuction Handle	Aluminum		0	